## **ALOE GEL SCAR FIXER**

Say goodbye to scars with Aloe Gel Scar Fixer, a potent aloe vera formula that fixes scars while soothing your skin, leaving you with a renewed sense of confidence.

Item	Ingredient	Wt%
1	Water	42.95
2	Triethanolamine	3.00
3	ALOE ACETYPOL POWDER (ACEMANNAN 15% TO 20%)	1.00
4	Chitosan	1.50
5	Onion extract	1.00
6	Dimethicone	40.00
7	Cetyl alcohol	3.00
8	Cellulose gum	2.00
9	Disodium EDTA	1.00
10	Potassium sorbate	0.50
11	Jojoba oil	2.00
12	Citric Acid	0.05
13	Propylene glycol	1.00
14	Methyl paraben	1.00



## Procedure

- 1. In a vessel heat the water between 50 60 °C, add the triethanolamine and aloe acetypol powder. Mix until is free of lumps.
- 2. Add the chitosan and onion extract to the mix keeping a temperature between 40-50 °C
- 3. In aseparate vessel heat the dimethicone around  $60 70^{\circ}$  with the ingredients from 7 to 12 with constant agitation at medium speed. 4. Mix both vessels in one batch, mix until is free of lumps.
- 5. Let cool at room temperature and add the conservatives (Propylene glycol and methyl paraben) when it's around 35 40 °C.

## \*DISCLAIMER

The information contained herein is being furnished for informational purposes only, upon the express condition that the User makes its own assessment of the appropriate use of such information. While the information contained herein is believed to be reliable, no representations, guarantees or warranties of any kind are made as to its accuracy, suitability for a particular application or the results to be obtained herefrom. Agromayal Botanica SA de CV. ("AMB") cannot guarantee how any products associated with this information will perform in combination with other substances or in User's process. Due to the variations in methods, conditions and equipment used commercially in processing the materials, on warranties or guarantees are made as to the suitability of the information/products for the applications disclosed. AMB-shall not be liable and the User assumes all risk and responsibility for, any use or handling of any material beyond AMB's direct control. AMB WELLLNESS MAKES NO WARRANTIES, EXPRESS OR INVIDENCE. It is the User's sole responsibility to determine if there are any issues relating to patent infringement of any component or combination of components relating to the supplied information. Nothing contained herein is be considered as permission, recommendation, nor as an inducement to practice any patented invention without the permission of the patent owner. This suggested for extensive testing for performance, efficacy or safety. BEFORE COMMERCIALIZATION, YOU SHOULD THOROUGHLY TEST. THE FORMULATION OF ANY VARIATION OF IT TO DETERMINE ITS PERFORMANCE, EFFICACY AND SAFETY. IT IS YOUR RESPONSIBILITY TO OBTAIN ANY NECESSARY GOVERNMENT CLEARANCE, LICENSE OR







