## ALOE SILK SCAR GEL

Experience the silky smoothness of Aloe Silk Scar Gel, a luxurious blend of aloe vera and botanical extracts that glides on effortlessly, providing effective scar treatment for visibly flawless skin.

	Item	Ingredient	Wt%	
	1	Water	39.55	
	2	Triethanolamine	3.00	
	3	<b>ALOE ACETYPOL POWDER (ACEMANNAN 15% TO 20%</b>	<u>6)</u> 1.00	
	4	Chitosan	1.50	
	5	Onion extract	1.00	
	6	Dimethicone	40.00	
	7	Cetyl alcohol	3.00	
	8	Oat lipid	4.00	
	9	Caprylic acid	0.70	
	10	Dissodium EDTA	0.25	
	12	Butylene Glycol	2.00	
	13	Mineral Oil	4.00	



## Procedure

- 1. In a vessel heat the water between 50 60 °C, add the triethanolamine and aloe acetypol powder. Mix until is free of lumps.
- 2. Add the chitosan and onion extract to the mix keeping a temperature between 40-50 °C
- 3. In separate vessel heat the dimethicone around  $60 70^{\circ}$  with the ingredients from 7 to 10 with constant agitation at medium speed. 4. Mix both vessels in one batch, mix until is free of lumps.
- 5. Let cool at room temperature and add the conservatives (Propylene glycol and methyl paraben) when it's around 35 40 °C.

## \*DISCLAIMER

The information contained herein is being furnished for informational purposes only, upon the express condition that the User makes its own assessment of the appropriate use of such information. While the information contained herein is believed to be reliable, no representations, guarantees or warranties of any kind are made as to its accuracy, suitability for a particular application or the results to be obtained herefrom. Agromayal Botanica SA de CV. ("AMB") cannot guarantee how any products associated with this information will perform in combination with other substances or in User's process. Due to the variations in methods, conditions and equipment used commercially in processing these materials, no warranties or guarantees are made as to the suitability for the information/products for the applications disclosed. AMB shall not be liable and the User assumes all risk and responsibility for, any use or handling of any material beyond AMB's direct control. AMB WELLLNESS MAKES NO WARRANTIES, EXPRESS OR INTEREST OF The IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. It is the User's sole responsibility to determine if there are any issues relating to patent infringement of any component or combination of components relating to the supplied information. Nothing contained herein is to be considered as permission, recommendation, nor as an inducement to practice any patented invention without the patent owner. This suggested formation is based are reliable, but it has not been subjected to extensive testing for performance, efficacy or safety. BEFORE COMMERCIALIZATION, YOU SHOULD THOROUGHLY TEST THE FORMULATION OF IT TO DETERMINE ITS PERFORMANCE, EFFICACY AND SAFETY. IT IS YOUR RESPONSIBILITY TO OBTAIN ANY NECESSARY GOVERNMENT CLEARANCE, LICENSE OR RESPONSIBILITY TO OBTAIN ANY NECESSARY GOVERNMENT CLEARANCE, LICENSE OR RESPONSIBILITY TO OBTAIN ANY NECESSARY GOVERNMENT CLEARANCE, LICENSE OR RESPONSIBILITY TO OBTAIN ANY NECESSARY GOVERNMENT CLEARANCE, LICENSE OR







